



**ASA DENTAL**  
MAKE PEOPLE SMILE

# DESIGN INSPIRED BY NATURE, CREATED TO EXCEL

AsaWave is the perfect meeting of art and science, combining the elegance of natural design with the precision of modern engineering. Born from the intuition of a visionary entrepreneur and the innovation of expert designers, AsaWave redefines the user experience of dental instruments.

Every detail is the result of careful research and development work, which combines creativity and functionality to offer a tool unique in its kind. Behind every curve there is an in-depth study that guarantees not only a

refined design but also unparalleled ease of use. Inspired by biomimicry, it studies the structural and functional principles present in nature, translating them into innovative solutions that enhance aesthetics and performance.

The fluid and harmonious shapes evoke the perfection of waves and natural contours, while functionality is optimized to ensure maximum practicality and efficiency. AsaWave is not just an instrument, but a natural extension of the professional, capable of transforming every procedure into a gesture of pure perfection.

## MORE THAN A PRODUCT: AN EXPERIENCE.

Every curve of AsaWave is designed to naturally adapt to the hand, becoming a fluid extension of your movements. The organic ergonomics optimize every gesture with naturalness, while the continuous surfaces repel dirt, facilitating immediate cleaning and preserving the purity of the design.



### MOUTH MIRROR HANDLE



RC2101

### SINGLE-ENDED EXPLORERS



6  
RC0700-6



9  
RC0700-9



23  
RC0700-23

### INSTRUMENT KITS FOR FIRST EXAMINATIONS



#### SRC0600-6

|          |                                           |   |
|----------|-------------------------------------------|---|
| RC0700-6 | explorer with AsaWave handle              | 6 |
| 0601-2   | London College cotton and dressing pliers | 2 |
| RC2101   | handle for AsaWave mouth mirror           |   |
| 2200-4   | plane mouth mirror                        | 4 |

#### SRC0600-9

|          |                                           |   |
|----------|-------------------------------------------|---|
| RC0700-9 | explorer with AsaWave handle              | 9 |
| 0601-2   | London College cotton and dressing pliers | 2 |
| RC2101   | handle for AsaWave mouth mirror           |   |
| 2200-4   | plane mouth mirror                        | 4 |

#### SRC0600-23

|           |                                           |    |
|-----------|-------------------------------------------|----|
| RC0700-23 | explorer with AsaWave handle              | 23 |
| 0601-2    | London College cotton and dressing pliers | 2  |
| RC2101    | handle for AsaWave mouth mirror           |    |
| 2200-4    | plane mouth mirror                        | 4  |

DOUBLE-ENDED  
EXPLORER



23/17  
RC0704-23D1

PERIODONTAL POCKET  
PROBES



CP12  
RC0702-12



Williams  
RC0703-14




Goldman-Fox  
Williams  
RC0707-GFW


ROOT CANAL  
PLUGGERS



1 Machtou  
RC0803-1



2 Machtou  
RC0803-2



3 Machtou  
RC0803-3

ROOT CANAL  
PLUGGERS




4 Machtou  
RC0803-4


EXCAVATORS



1  
RC1700-1



2  
RC1700-2




3  
RC1700-3


HEIDEMANN SEPARATING  
SPATULAS



1  
RC1102-1



2  
RC1102-2




3  
RC1102-3

AMALGAM FILLING INSTRUMENTS

PLASTIC FILLING  
INSTRUMENTS

FILLING INSTRUMENTS  
Burnishers

COMPOSITE INSTRUMENTS





GRACEY USE CHART

- Anterior sextants of teeth

GRA 1/2

Incisors and Canines

GRA 3/4

Incisors and Canines

GRA 5/6

Anteriors and Premolars
- Buccal and lingual surfaces of posterior teeth

GRA 7/8

Premolars and Molars

GRA 9/10

Premolars and Molars
- Distal surfaces of posterior teeth

GRA 11/12

Premolars and Molars

GRA 15/16


Premolars and Molars
- Mesial surfaces of posterior teeth

GRA 13/14

Premolars and Molars

GRA 17/18

Molars




GRACEY CURETTES




Set of 7  
SRC1807




1/2 Gracey  
RC1807-02




3/4 Gracey  
RC1807-04




5/6 Gracey  
RC1807-06




7/8 Gracey  
RC1807-08




9/10 Gracey  
RC1807-10




11/12 Gracey  
RC1807-12



13/14 Gracey  
RC1807-14



15/16 Gracey  
RC1807-16



17/18 Gracey  
RC1807-18

GRACEY CURETTES


SHARPEN FREE

A new type of coating: **DLC Diamond-Like Carbon**.  
DLC is known for its **high degree of hardness**, ensuring greater wear resistance. DLC has a **low coefficient of friction**, reducing residue adhesion and heat generation during instrument use.  
The good **corrosion resistance** protects dental instruments from oxidation and wear caused by chemicals in the mouth.  
DLC-coated dental instruments can be **easily sterilized** with standard methods, such as autoclaving, without compromising the coating properties.  
Our DLC-coated instruments offer **extended durability, consistent precision, and enhanced safety**. This means longer-lasting tools that save you time and reduce costs, providing a tangible investment in your practice. You'll deliver **high-quality treatments** with optimal patient **comfort and safety**.



CURETTES AND SCALERS

SHARPEN FREE




CURETTES




13/14  
Columbia  
RC1806-14




13/14 Mc Call  
RC1825-14




6 Crane-Kaplan  
RC1829-CK6



6S mini  
Crane-Kaplan  
RC1829-CK6S



204S  
RC1827-204S




13S/14S Mc Call  
RC1825-14S


SCALERS




M23  
universal scaler  
RC1900-M23




U15/33  
RC1902-1



U15/30  
Towner-Jacquette  
RC1902-2




H5/33  
RC1904-H5



H6/H7  
RC1904-H6



BONE CURETTES



85 Lucas  
RC2001-85

86 Lucas  
RC2001-86

87 Lucas  
RC2001-87

88 Lucas  
RC2001-88







#### HEAD OFFICE

##### **Asa Dental S.p.A.**

Via Valenzana di Sopra, 60 - Bozzano  
55054 Massarosa (Lu) - Italy  
Tel. +39 0584 938305 - 93363 - 938306  
Fax +39 0584 937167  
asadental@asadental.it

#### LOGISTICS CENTER

Via Sarzanese Valdera, 1331 - Bozzano  
55054 Massarosa (Lu) - Italy

#### PRODUCTION PLANTS

Via delle Piagge, 1/A - Martia  
55014 Capannori (Lu) - Italy

Via Fierla, 13 - Z.I.  
33085 Maniago (Pn) - Italy

#### BRANCH OFFICES

##### **Asa Dental USA Corp.**

105 Maxess Road, Suite 124  
Melville, New York 11747 - USA

##### **Asa Dental CHINA CO. LTD.**

Room 315, Huaying Building,  
Zhong Xin Avenue  
Tianjin Airport Economic Area  
Tianjin - China

##### **Asa Dental S.p.A. (Branch)**

Shahrazad Afdal Bahati Bldg.  
Parcel 345-826  
Bur Dubai  
PO Box 414299  
Dubai - UAE

**www.asadental.com**

